

 **cuattro DR^{HD}**



A whole new category of Digital Radiography...


Cuatro DR HD (High Definition) Digital Radiography is now available exclusively from scil. Cuattro DR HD delivers detail, resolution, and contrast image quality that is visibly better, much better, than Standard Definition systems.

cuattro DR^{HD}



Optional Table Specifications:

- Generator: 32kW, 400mA, 125kvp
- 4-way Float 60" x 30" Table Top
- Software Controlled Generator
- Medical Grade LED Collimator



HD Minimum Standards:

- Direct Deposition Cesium **Cs**
- 100 micron pixel size
- HD pixel density. 18 Million (17x17)
- HD pixel density. 15 Million (14x17)
- Full Size Image Area (14x17 - 17x17)
- Resolution of at least 5 lp/mm
- 16 Bit A/D 16 Bit Grayscale Native
- 16 Bit Multi-Frequency Tuning

Cuattro DR HD Added Benefits:

- Automatic Exposure Detection
- Cassette Sized Form Factor
- Embedded - 64 bit OS
- Multi-Touch User Interface
- Cloud Based PACS Connectivity
- Cloud Archival
- Software Upgrades and Updates
- Canadian Service and Support

A smarter direction in Digital Radiography. Better, much better.

www.heskavet.ca • toll free: +1.866.382.6937 • imaging the possibilities

© 2024 Heska Corporation. All rights reserved. Antech, the Antech logo, and all other trademarks used herein are the registered trademarks of Antech Diagnostics, Inc. or its affiliates. HESKA is a registered trademark and Cuattro DR HD is a trademark of Heska Corporation in the United States and in other countries. Specifications are subject to change without notice and may vary based on operating conditions. Pricing dependent on final configuration. CA20LT0602_005



HESKA
AN ANTECH® COMPANY