



Vet/IV 2.2

Infusion Pump

Accurate, portable infusion therapy.

Find out why thousands of customers love the HESKA® VET/IV™ 2.2 Infusion Pump. Weighing only 2.6 pounds, the Vet/IV Pump supports three pressure settings and can be programmed for very small to very large volumes with rate and time or rate and volume to be infused. With multiple programming options, intuitive 6 button interface, and communication alarms for any issue, Heska's Vet/IV 2.2 Infusion Pump is a reliable necessity in the veterinary practice.

Lightweight and convenient

Battery powered

Trusted by thousands

Intuitive operation

Easy Operation

Use for fluid replacement or to administer total parenteral nutrition (TPN) even in very small doses.

Features:

Lightweight and Portable:

- Weighs only 2.6 pounds and fits in the palm of your hand.
- Operate off rechargeable battery or AC power.

Intuitive Operation:

- Set all functions using just six buttons.

Convenient:

- Hang from a treatment cage or IV pole.

Vet/IV



World
Class
Support

Heska strives to ensure that every customer experience is world class. Each Heska instrument is supported by highly trained, dedicated technical specialists.

Reliability

Never be unsure of your Vet/IV Pump status with Communication Alarm settings that alert to issues with low battery, end battery, air in line, high pressure, fill set, replace set and more.



Flexible Programming Options

Programmable for very small to very large volumes with rate and time or rate and volume to be infused.

Multiple calibration options with fluid line sets available for:

- Abbott® Series A (same calibration as Hospira® IV sets)
- Abbott Series B (same calibration as Baxter® IV sets)
- Baxter ACT 5435s
- Hospira List No. 11545–58

Best-In-Class Technology

Micro and Macro Flow Rate:

- Flow rate ranges from 0.1–999 ml/hour.

Multi-purpose:

- Use for fluid replacement even in very small doses.

Multiple Programming Options:

- Program up to 10 drip rates and times.
- Program rate and time or rate and volume.

Multiple Pressure Settings:

- Select from three pressure settings for different patient conditions.

Piggyback Capable:

- Option to control the flow of fluid from two fluid bags.

www.heskavet.ca • 866.382.6937

©2024 Heska Corporation. All Rights Reserved. Antech, the Antech logo, and all other trademarks used herein are the registered trademarks of Antech Diagnostics, Inc. or its affiliates. HESKA is a registered trademark and Vet/IV is a trademark of Heska Corporation in the U.S. and other countries. CA20LT0709



HESKA®

AN ANTECH® COMPANY