

NEW



Therapy Lasers

For acute or chronic, superficial or deep conditions, laser therapy can help. By accelerating the healing process with increased circulation and enhanced tissue repair, laser therapy decreases inflammation and pain. In the management of chronic and acute pain, Class IV lasers created by Summus Medical Laser® are a preferred choice for practitioners worldwide.

Condition Specific Protocols

Portable & Battery Operated

Wavelength Selectability

H1 Laser

With a lightweight yet sturdy design, the H1 delivers powerful results and simple operation. It's built on an advanced Android platform for easy and quick selection of an expert treatment protocol that will get your patients healing better, faster.

- 165 treatment protocols
- Create and save custom protocols
- Portable and sturdy
- Lightweight yet durable construction

3 wavelengths

- 915 nm enhances oxygen delivery
- 810 nm increases ATP production
- 650 nm accelerates surface healing



P4 Laser

Not only the most advanced Class IV therapy laser, but also built on a technology platform that redefines the therapy laser industry. No other therapeutic laser can deliver usage metrics and tools on the laser interface to help your practice succeed with laser therapy.

- Faster treatment times
- 1,300+ treatment protocols
- Programmable custom protocols

4 wavelengths

- 980 nm improves circulation
- 915 nm enhances oxygen delivery
- 810 nm increases ATP production
- 650 nm accelerates surface healing
- The right combination of power, wavelengths and pulsing delivers an optimal therapeutic dose in less time

5 hand piece lenses

- Maximum treatment versatility with our 50 mm broad beam (open lens and curved contact lens), 25 mm zoom (open lens), 25 mm (curved contact lens) and ENT lenses.



Laser Therapy Aids:

- Osteoarthritis
- Joint pain
- Tendinopathies
- Edema and congestion
- Ligament sprains
- Muscle strains
- Puncture wounds
- Post-traumatic injury
- Post-surgical pain
- Neck and back pain
- Hip dysplasia
- Burns
- Chronic wounds
- Rehabilitation
- Post-orthopedic surgical recovery

LT

Stay Connected

Summus Connect Software

- Patient modality worklist
- Add patient notes
- Reporting to maximize ROI
- Generate invoices
- Generate appointments and see them on your laser

Only available on the Horizon Laser System, Summus Connect Software seamlessly integrates your new laser with any practice management software that accepts DICOM images! Imagine being able to track all of your laser patients in one place and to generate reports in a totally paperless environment

The days of missed charges and duplicate entries are over. Have your laser therapy workflow, and your practice, running smoother than ever, with Connect.

Horizon Laser

The first of its kind complete laser system that allows you to easily and affordably integrate laser therapy into your practice. We partner with you to provide consistent, precise, reproducible results for your patients.

Independent wavelength control

- Each wavelength has a different effect on the body. This allows a condition to receive the optimal therapy.

Automatic optic identification

- The laser automatically recognizes which optic is being used. This is important from a safety and user experience standpoint.
- Progressive treatments
- A patient can be set-up for a treatment plan. The protocol will vary as the patient progresses through the plan.

4 wavelengths

- 980 nm improves circulation
- 915 nm enhances oxygen delivery
- 810 nm increases ATP production
- 650 nm accelerates surface healing

Phased protocols

- Each protocol is built with multiple phases of wavelengths and pulsing.

Visual guidance

- Motion graphics displayed on-screen will guide you through treatments with confidence.

