

AESCULAP® Acculan 4 Power System

A powerful and simple to use system which combines reliability, longevity and ergonomics to meet all your veterinary surgical needs. The rounded organic shape of the grip area provides ergonomic comfort, and also serves as a functional housing for the battery.

One battery for two systems

Lightweight and ergonomic



Battery integrated handle, reducing bulk

AESCULAP® Acculan 4 Drill

The German engineered Acculan 4 lithium ion battery operated drill system offers performance, the highest standards of reliability, durability and ergonomics. Its lightweight design offers versatility, and the fast attachment system allows for a wide range of surgical procedures.



Available Attachments



AESCULAP® Acculan 4 TPLO Saw

The Acculan 4Ti TPLO Saw is the next generation of battery operated saw system made for veterinarians. This saw is designed to deliver an accurate cut with maximum power.

The ergonomic design offers easy handling with surgical precision.

The novel motor and control concept offers the kind of performance previously found only in the Elan and pneumatic systems. This saw is only for TPLO procedures.

Features

- Titanium housing
- Ergonomic design for a stable cut
- Patented Aesculap motor
- Contactless mechanism
- Variable speed
- Continuous start
- Powerful torque
- Standard triangular shaft blades
- Battery sufficient for more than 10 TPLO's
- Non-sterile battery

One Battery for Two Systems



LI-ION Technology

- More power than previous NiMH batteries
- Longer life cycle
- No memory effect
- Reliable



www.heskavet.ca 866.382.6937

©2023 Heska Corporation. All Rights Reserved. Aesculap is a registered trademark of Aesculap is a registered trademark of Aesculap AG. HESKA is a registered trademark of Heska in the U.S. and other countries. CA23LT0604

