

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0524A



2024-07-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0524A-1	LOT	ME0524A-2	LOT	ME0524A-3
element HT5+	Ret# x10 ³ /μL	84.0	± 64.4	230.2	± 99.0	192.9	± 99.0
QC Mode	Ret%	1.50	± 1.12	4.93	± 2.00	8.65	± 3.50
	IRF %	13.5	± 15.0	9.8	± 13.0	7.5	± 25.0
	LFR %	86.5	± 15.0	90.2	± 15.0	92.5	± 25.0
	MFR %	12.3	± 12.0	7.7	± 10.0	4.3	± 22.0
	HFR %	1.2	± 10.0	2.1	± 10.0	3.2	± 10.0
	RHE pg	26.3	± 12.0	25.9	± 12.0	21.3	± 12.0
	IPF %	7.4	± 15.0	7.7	± 15.0	7.7	± 15.0

Before using, refer to the instruction sheet for mixing directions.